

Development of High-Glucosinolate-Retaining Lactic-Acid-Bacteria-Co-Fermented Cabbage Products

Chen-Che Hsieh 1, Yuan-Heng Liu 2, Shin-Ping Lin 3,4,5,6 ORCID, Shella Permatasari Santoso 7 ORCID, Kaemwich Jantama 8 ORCID, Tsung-Yu Tsai 9, Chang-Wei Hsieh 10, *ORCID and Kuan-Chen Cheng

Abstract

This study enhanced cabbage nutrition and functionality via bioreactor-based lactic acid fermentation, focusing on glucosinolate retention. Using *Lactiplantibacillus plantarum*, *Lactobacillus acidophilus*, and *Bifidobacterium longum*, fermentation (35 °C, 24 h, 5 rpm) significantly improved antioxidant activity (16.32% increase in DPPH scavenging). Glucosinolate retention reached 82.02%, and remained at 66.49% after 14 days of storage at 4 °C. Results support bioreactor fermentation as a promising method for boosting cabbage's nutritional and commercial value.



Winpact Model: FS-V-SA05P

Introduction

Cabbage (*Brassica oleracea var. capitata*) is a widely cultivated, low-cost vegetable rich in glucosinolates, known for antioxidant and anticancer effects. However, these compounds degrade during processing via myrosinase activity. Lactic acid fermentation lowers pH, inhibiting this enzyme and preserving glucosinolates. Solid-state fermentation is cost-effective and environmentally friendly. This study used co-fermentation with *Lactiplantibacillus plantarum*, *Lactobacillus acidophilus*, and *Bifidobacterium longum*, followed by scale-up in a 5 L bioreactor to enhance glucosinolate retention, storage stability, and probiotic value.

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai, China
Postal code: 200333
Cell: +86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2, Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
F/ +91-40-27001586

Materials and Methods

For the scaled-up production trial, a 5 L solid-state bioreactor (FS-V-SA05P, Major Science, Taoyuan, Taiwan) was used. The parameters were as follows: the fermentation temperature was set at 35 °C, the fermentation time was 24 h, the stirring speed was 5 rpm, and the inocula of lactic acid bacteria (10^7 CFU/mL) were 0.3% and 3.0% (w/w) of the weight of the salted dehydrated cabbage.

Results

This study implemented a 5 L bioreactor-based co-fermentation of cabbage using *B. longum* BCRC 14634, *L. acidophilus* BCRC 14079, and *L. plantarum* BCRC 11697 to retain antioxidant properties and glucosinolates. Fermentation increased free phenolic and total flavonoid contents by 41.13% and 24.44%, respectively. The lowered pH inactivated myrosinase, achieving 82.02% glucosinolate retention. Under 35 °C and 5 rpm, even with reduced inoculum size, retention remained high. After 14 days at 4 °C, glucosinolate content and LAB count declined by 23.85% and 4.40 log CFU/g. The method showed superior retention over traditional processing.

References

Development of High-Glucosinolate-Retaining Lactic-Acid-Bacteria-Co-Fermented Cabbage Products

<https://www.mdpi.com/2311-5637/10/12/635>

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai, China
Postal code: 200333
Cell:+86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586



www.majorsci.com
info@majorsci.com