

NEW

ODM-MINI-500-CP

# Power Supply

MS major science  
www.majorsci.com  
info@majorsci.com



## Product Introduction

The Mini-500CP Volt Power Supply is both compact and versatile, making it the perfect choice for those with limited funds or space for lab equipment.

Fixed voltage or current can be selected according to the experimental needs. The compact size and complete functions meet your basic experimental needs, such as horizontal or vertical electrophoresis.

## Features

- ◆ New product design
- ◆ Operates at constant voltage or current
- ◆ Voltage and current settings with a minimum adjustment scale of 1 volt or 1 milliamp.
- ◆ Output voltage: 500 volts
- ◆ Output current: 400 mA
- ◆ Output power: 120 watts
- ◆ 4 output terminals
- ◆ Timed warning function
- ◆ Safety device: no load detection
- ◆ Wrapped plugs and sockets
- ◆ Excellent price/performance ratio
- ◆ With horizontal or vertical electrophoresis

## Specification

Output Voltage / Inc.	Step1: 10 - 500V / 1V Step2: 10 - 500V / 1V
Output Current / Inc.	Step1: 1 - 400mA / 1mA Step2: 1 - 400mA / 1mA
Max. Watt	120W
Output Type	Constant voltage or constant current
Control	Microprocessor controller
Terminal Pairs	4 pairs
Timer	999 (min) with alarm/ Continuous (set=0)
Safety Device	No load detection; shrouded plugs and sockets
Operation Temperature	4°C-40°C
Dimension	156 x 212 x 84mm (W x L x H)
Material	Polycarbonate housing and aluminum base plates
Weight	Approx. 1.1 kg
Rated Voltage	100 - 240 V~, 50-60Hz, 2.5A
Program	2-STEP Program

## Ordering Info

Cat. No.	Description
ODM-MINI-500-CP	ODM Mini Pro 500V Power Supply, CP version

\* All images are for reference only, actual products might differ from the pictures above.  
\* Technical specifications subject to change without notice.

### Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,  
Taoyuan City 330041, Taiwan  
T/+886-3-3762878  
F/+886-3-3761310

### US Office

19959 Sea Gull Way  
Saratoga, CA 95070  
U.S.A.  
T/ +1-408-366-9866  
F/ +1-408-446-1107

### Shanghai Office

Room 612, International business exhibition center,  
9300 Hunan Road, Pudong, Shanghai, China  
National toll-free No.:400-823-9177  
T/ +86-21-50795277  
F/ +86-21-50795277

### India Office

D.No.7-143, 2nd Floor,  
St.No.2,Nagendra Nagar,  
Habsiguda, Hyderabad-500007.  
India  
T/ +91-40-27001515  
T/ +91-40-27001586