



Innovative Life Sciences Tools



major science Tip

- **Universal Compatibility**

Standard Size Design: Compatible with a wide range of pipette brands commonly found on the market.

- **Operational Accuracy**

Ultra-Smooth Inner Wall Design: Effectively reduces residual liquid and enhances pipetting precision.

Optimized Orifice Design: Ensures smooth liquid flow without clogging.

- **Material and Safety**

Manufactured with Virgin PP Material: Complies with USP Class VI standards; excellent biocompatibility.

Filter Made of Ultra-High Molecular Weight Polyethylene (UHMW-PE): Ultra-pure, contaminant-free, and highly hydrophobic.

- **Contamination Prevention**

Hydrophobic Filter Design: Forms a barrier against aerosols to effectively prevent cross-contamination between samples and pipettes.

Clean Production Process: Free from DNase, RNase, and pyrogens.

- **Thermal Resistance & Molding Technology**

Wide Temperature Resistance: -80°C to 121°C; maintains shape after high-temperature sterilization.

Ultra-Fine Polished Molding Without Release Agent: Reduces chemical residues and ensures stable product quality.

Cat No.	Product Description	Packaging Method
MS-PW-T10SA1	10µl Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-T10A1	10µl Tip	1000pcs/bag, 20bag/carton
MS-PW-TE10SA1	10µl Extended Length Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-TE10A1	10µl Extended Length Tip	500pcs/bag, 2bag/box, 10box/carton
MS-PW-TFE10SA1	10µl Extended Length Filter Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-TF10SA1	10µl Filter Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-T200SA1	200µl Tip,	96pcs/box, 50box/carton, Sterilizable
MS-PW-T200A1	200µl Tip	1000pcs/bag, 20bag/carton
MS-PW-TE200SA1	200µl Extended Length Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-TE200A1	200µl Extended Length Tip	1000pcs/bag, 20bag/carton
MS-PW-TF200SA1	200µl Filter Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-TFE200SA1	200µl Extended Length Filter Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-TE1250SA1	1250µl Extended Length Tip	96pcs/box, 50box/carton, Sterilizable
MS-PW-TE1250A1	1250µl Extended Length Tip	1000pcs/bag, 5bag/carton
MS-PW-TFE1250A1	1250µl Extended Length Filter Tip	96pcs/box, 50box/carton, Sterilizable

* All images are for reference only, actual products might differ from the pictures above.

* Technical specifications subject to change without notice.



US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai
Postal code: 200333
Cell: +86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2, Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586

