

NEW


Winpact CO₂ resistant orbital shaker with 12" colorful touch screen is developed to offer an integrated solution for every mixing requirement.


WP
Winpact


Winpact CO₂ resistant orbital shaker, WP-COS

- Corrosion-resistant construction protects internal components from high CO₂ and humidity.
- Compact design enables use inside a CO₂ incubator for versatile applications.
- Platform supports multiple accessories, including flask clamps and tube racks, to accommodate lab labware.
- Variable speed range of 1–300 rpm ensures smooth shaking.
- 12" touch screen for intuitive operation.

Cat. No.	WP-COS
Speed	0, 1 to 300 rpm (1 rpm step)
Motor	Brushless motor
Top Plate Material	316 L
Capacity	Max. 5.2 kg
Controller	12" color touch panel
Mode	Continuous / Timing / Programable
Voltage	AC 100-240V, 50/60Hz
Dimensions (W x L x H)	Controller : 320 x 102 x 246 mm Base : approx. 310 x 310 x 97 mm (including tray platform)
Weight	Controller : approx. 4.65 kg / Base : approx. 7.15 kg

Cat. No.	Product Description
WP-COS	CO ₂ resistant orbital shaker with 12" colorful touch screen
WP-12-TC-CONTROL	12" colorful touch screen
WP-COS-BASE	Shaker base unit of Tetrad Slow Speed Magnetic Stirrer with 12" colorful touch screen
WP-12-TC-CABLE	Ribbon Cable set for 12" colorful touch screen use

major science
Innovative Life Sciences Tools



* All images are for reference only, actual products might differ from the pictures above.
* Technical specifications subject to change without notice.

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Shanghai Office

Building #C1, 880,
JiangYang Nan Road,
Baoshan, Shanghai,
china 200439
T/ 13585666474

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

India Office

D.No.7-143, 2nd Floor,
St.No.2, Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586