

New

Tetrad Slow Speed Magnetic Stirrer WP-SMS

WP
Winpact



- ▶ 12" touch screen for user-friendly observation.
- ▶ Hermetically SUS316L housing is easy to clean and resistant to exposure to CO₂ environments.
- ▶ Smooth and even stirring at speeds as low as 1 rpm.
- ▶ Four independent stirring positions, each capable of handling up to 1 L of culture volume simultaneously.
- ▶ Each stirrer can be operated individually, allowing flexible and separate control across all positions.
- ▶ A Continuos / Timing / Programmable stirring mode is available.

Cat No.	WP-SMS
Speed	0, 1 to 300 rpm (1 rpm step)
Stirring Positions	4 positions, individual rotation control
Top Plate Material	316L
Capacity	Max. 4 x 1L 2 x 2L(diagonal position) 2 x 3L(diagonal position)
Controller	12" color touch panel
Mode	Continuous / Timing / Programable
Voltage	AC 100-240V, 50/60Hz
Dimensions (W x L x H)	Controller : approx. 320 x 102 x 246 mm Base : approx. 310 x 310 x 46 mm
Weight	Controller : approx. 4.65kg Base : approx. 4.4kg

Cat No.	Product Description
WP-SMS	Tetrad Slow Speed Magnetic Stirrer with 12" colorful touch screen
WP-12-TC-CONTROL	12" colorful touch screen
WP-SMS-BASE	Stirrer base unit of Tetrad Slow Speed Magnetic Stirrer with 12" colorful touch screen
WP-12-TC-CABLE	Ribbon Cable set for 12" colorful touch screen use



MS major science
www.majorsci.com
info@majorsci.com

4 positions,
individual rotation control

* All images are for reference only, actual products might differ from the pictures above.
* Technical specifications subject to change without notice.



US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai
Postal code: 200333
Cell:+86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586